

INSTRUMENTS TABLE:

Material of Construction	
Under Bench Modules	<p>Under bench Modules will be made of 20 SWG pre-coated Galvanised iron Sheet (TATA Steel - GI Sheet) and completely finished with epoxy powder coating.</p> <p>These modules will be provided with top drawer and bottom cupboard with removable inner partition and with lockable doors and aesthetic powder coated handles.</p> <p>The under-bench modules will be arranged alternatively. The base of the modules will be provided with skirting panels for aesthetic look.</p>
Knee hole Space	<p>Knee hole space will be provided as per the diagram.</p>
Colour	<p>Colour combination of IVORY and drawer fronts contrasting BLUE for aesthetic appeal by powder coating technique.</p>
Granite Slab	<p>Made of finely polished 18 mm \pm 2 mm thick black Colour Granite with rounded edges.</p>
Electrical points	<p>15 / 5 Amps Electrical Points will be provided as per the drawing.</p>

