

## WALL TABLE:

<b>Material of Construction</b>	
<b>Under Bench Modules</b>	<p>Under bench Modules will be made of 20 SWG pre-coated Galvanised iron Sheet (TATA Steel - GI Sheet) and completely finished with epoxy powder coating.</p> <p>These modules will be provided with top drawer and bottom cupboard with removable inner partition and with lockable doors and aesthetic powder coated handles.</p> <p>The under-bench modules will be arranged alternatively. The base of the modules will be provided with skirting panels for aesthetic look.</p>
<b>Knee hole Space</b>	<p>Knee hole space will be provided as per the diagram.</p>
<b>Colour</b>	<p>Colour combination of IVORY and drawer fronts contrasting BLUE for aesthetic appeal by powder coating technique.</p>
<b>Granite Slab</b>	<p>Made of finely polished 18 mm <math>\pm</math> 2 mm thick black Colour Granite with rounded edges.</p>
<b>Sink</b>	<p>PP Sink with swan neck tap and waste pipe up to lower portion of the table will be provided. - quantity as per the diagram</p>

